

Technique Specification

Physical Specifications

Dimensions	50*46*32 cm
Weight	22 kg
Casters	Front wheels with lock

Pipeline Supply

Gas configuration	O <sub>2</sub> , N <sub>2</sub> O
Pipeline input range	280 – 600 kPa

Modes of Ventilation

MANUAL	Manual
--------	--------

Ventilator Parameter

Working Type	Pneumatic driven
--------------	------------------

System Standard

Hypoxic guard system	N <sub>2</sub> O cut-off valve, O <sub>2</sub> concentration 25%
APL valve	< 6 kPa
O <sub>2</sub> Flush	25 ~ 75 L/min
Volume of CO <sub>2</sub> absorber	0.8 L (Optional for 1.6 L CO <sub>2</sub> absorber)
Flowmeter	4 tubes, Cascade 0.05 ~ 10 L/min O <sub>2</sub> , 0.05 ~ 10 L/min N <sub>2</sub> O
Vaporizer	1 unit, (Isoflurane / Enflurane / Sevoflurane / Halothane)

Optional

Anesthesia ventilator (CWH-1010, CWH-1020)
Trolley
AGSS



NANJING CHENWEI MEDICAL EQUIPMENT CO., LTD

- #93 Shengtai West Road, Jiangning, Nanjing, Jiangsu, China 211100
- +86 25 86989898-827/836
- sales@chenwei-med.com
- www.chenwei-med.com

E-catalog



CWM-101A  
VET Anesthesia Machine







## Safe, stable and fully functional suitable for various environments

### Feature

- Portable design, easy to operate and save space in operation theatre.
- Molding Production with a benefit of light weight and rust-proof.
- Integrated breathing circuit design ensures simple operation and clear equipment. 0.8 L CO<sub>2</sub> absorber also equipped on breathing circuit.
- APL valve design can set the pressure limit under manual control.
- 4 tubes flowmeter for O<sub>2</sub> and N<sub>2</sub>O, support low flow anesthesia.
- Self-designed high precision Vaporizer with flow, pressure, temperature compensation, self-locking structure. High precision of anesthesia agent control.



### VET Anesthesia System Connected To Veterinary 1-6# Mask Set

- Combination Instruction:
- Combination 1: 1#, 2# combined, adapt to round face animal, such as Cat and new born animal.
- Combination 2: 3#, 4#, 5#, 6# combined, adapt to long face animal, such as Dog, Pig, Rabbit, Bird.
- Combination 3: 1#, 2#, 6# combined, adapt to small and big Mouse.
- Combination 4: 1#~6# combined, full size, wide application.

- Combination 1:   ...
- Combination 2:     ...
- Combination 3:   ...



### Available To Connect Anesthesia Ventilator (CWH-1010, 1020)

- Multi Ventilation modes applicable to all kinds of animals.
- Suitable for all kinds of animals below 200 kg.
- Portable design ventilator can directly connect to anesthesia machine.

