



Professional manufacturer for Anesthesia & Ventilator Since 1992

NANJING CHENWEI MEDICAL EQUIPMENT CO., LTD

- #93 Shengtai West Road, Jiangning, Nanjing, Jiangsu, China 211100
- © +86 25 86989898-827/836
- sales@chenwei-med.com
- www.chenwei-med.com





Integrated Breathing Circuit

- Integrated breathing circuit with APL valve provides excellent tightness and safety.
- The circuit components resistant to high temperature at 134°C and high pressure disinfection.
- Replacing drive bellow is not required when operating adult / pediatric anesthesia.
- 2L volume sodium lime tank, single tank structure is more convenient to disassemble.





High Precision Vaporizer

- Vaporizer can support low flow anesthesia (as low as 0.2 L/min) and high flow anesthesia (up 15 L/min).
- Accurate output of concentrations under varying conditions of pressure, flow and temperature.
- Support Selectatec mounting type and cagemount compatible.
- Agent type: Isoflurane, Enflurane, Halothane, Sevoflurane.
- Maintenance-free.
- Fill-in type: pouring.

Yoke system (Optional)

Up to 3 back up cylinders with A-type 11 L size.







Features



- 10.1 inch Large screen.
- Intuitive English operating interface.
- Comprehensive monitoring capabilities include ventilation parameters monitoring, preset parameters, alarm indication and breathing waveform monitoring.
- Ventilation modes: VCV, SIMV, A/C, MANUAL, STAND-BY.

- Max. 2 high-precision vaporizers.
- FlowMeter: 2 tubes, O₂: 0.05 15 L/min;
 2 tubes, N₂O: 0.05 15 L/min;
 1 tube, Air: 0.1 15 L/min (Optional).





ACGO (Optional)

- International standard design.(diameter: 22mm).
- Convenient for patients'oxygen inhalation and postoperative recovery
- To connect with T circuit or bain circuit for open-operation.